CeraNews

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Brazil

Daniel Stainer
dstainer@terra.com.br
0055 / 48 / 99 78 49 42

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1. CeraNews 2/2008: 2
2. Kelley SS, Lachiewicz PF, Hickman JN, Paterno SM (1998): Relationship of femoral head and acetabular size to the prevalence of dislocation. *Clin Orthop* 355: 163-170
3. Beaver RJ, McCormick C, KopA (2002): Fractures in misaligned ceramic acetabular liners. In: Garino JP, Willmann G (Hrsg.). *Bioceramics in Joint Arthroplasty*. Thieme Verlag Stuttgart, NewYork: 65-71
4. Zagra L, Giacometti Ceroni R, Corbella M (2004): THA Ceramic-Ceramic Coupling: The Evaluation of the Dislocation Rate with Bigger Heads. In: Lazennec JY, Dietrich M (Hrsg.). *Bioceramics in Joint Arthroplasty*. Steinkopff Verlag, Darmstadt: 163-169

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1. Baek S-H, Kim S-Y, Cementless Total Hip Arthroplasty with Alumina Bearings in Patients, Younger Than Fifty with Femoral Head Osteonecrosis, *J Bone Joint Surg Am.* 2008;90:1314-1320
2. Tözün I, Erdem M, Dikici F, Alumina ceramic-on-ceramic total hip arthroplasty: 5-6 year follow-up study, *Poster EFORT 2008*
3. Hutten D, Lestrat V, Loriaut P, The alumina-alumina bearing for total hip arthroplasty in patients aged less than 60 years: results with eight years mean follow-up, *Poster EFORT 2008*
4. Cruz Pardos A, Garcia Rey E, Garcia Cimbrello E, Ortega Chamorro J, Alumina-on-alumina THR in patients with juvenile idiopathic arthritis: a five-year follow-up prospective study, *Poster EFORT 2008*
5. Garcia Rey E, Cruz Pardos A, Garcia Cimbrello E, Ortega Chamorro J, Alumina-on-alumina total hip arthroplasty in patients younger than 40 years. A prospective multicenter study, *Poster EFORT 2008*
6. O’Gorman A, Waters T, Walter W, Zicat B, Walter W, ABG II Ceramic-on-Ceramic Total Hip Arthroplasty – 5 Year Minimum Follow up, *Poster EFORT 2008*
7. Zagra L, Giacometti Ceroni R, Ceramic-on-ceramic coupling with 36 mm heads, *Poster EFORT 2008*
8. Zhinian Wan, Myriam Boutary, Lawrence D. Dorr, The Influence of Acetabular Component Position on Wear in Total Hip Arthroplasty, *Journal of Arthroplasty Vol. 23 / 1*
9. Needham, J, Burns T, and Gerlinger T, Catastrophic Failure of Ceramic Polyethylene Bearing Total Hip Arthroplasty, *Journal of Arthroplasty Vol. 23 No. 4/2008*

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- Y. Catonné, J.-Y. Lazennec, A. Nogier, E. Fourniols and B. Masson, *Revision strategy after metal on metal THR failure: Conversion to ceramic on ceramic*, in: Benazzo F, Falez F, Dietrich M (eds), *Bioceramics and Alternative Bearings in Joint Arthroplasty*, 11th BIOLOX® Symposium Proceedings, Steinkopff Verlag Darmstadt: S. 141–144, 2006
- T. Güttler, *Experience with BIOLOX® OPTION revision heads*, in: Benazzo F, Falez F, Dietrich M (eds), *Bioceramics and Alternative Bearings in Joint Arthroplasty*, 11th BIOLOX® Symposium Proceedings, Steinkopff Verlag Darmstadt: S. 149–154, 2006
- K. Knahr and M. Pospischill, *Strategies for Head and Inlay Exchange in Revision Hip Arthroplasty*, in: Chang JD, Billau K (eds), *Bioceramics and Alternative Bearings in Joint Arthroplasty*, 12th BIOLOX® Symposium Proceedings, Steinkopff Verlag Darmstadt: S. 275–280
- G. Pignatti, C. Stagni, D. Dallari, V. Bochicchio, A. Raimondi and A. Giunti, *Modular neck and ceramic on ceramic coupling in revision total hip arthroplasty*, in: Benazzo F, Falez F, Dietrich M (eds), *Bioceramics and Alternative Bearings in Joint Arthroplasty*, 11th BIOLOX® Symposium Proceedings, Steinkopff Verlag Darmstadt: S. 145–148, 2006

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1. Dreinhöfer KE, Dieppe P, Stürmer T, Gröber-Gratz D, Flören M, Günther KP, Puhl W, Brenner H (2006) Indications for total hip replacement: comparison of assessments of orthopaedic surgeons and referring physicians. *Ann Rheum Dis* 65:1346-1350
2. Hawker GA, Wright JG, Coyte PC, Williams JI, Harvey B, Glazier R, Wilkins A, Badley EM (2001) Determining the need for hip and knee arthroplasty: the role of clinical severity and patients' preferences. *Med Care* 39:206-216
3. Hudak PL, Clark JP, Hawker GA, Coyte PC, Mahomed NN, Kreder HJ, Wright JG (2002) "You're perfect for the procedure! Why don't you want it?" Elderly arthritis patients' unwillingness to consider total joint arthroplasty surgery: a qualitative study. *Med Decis Making* 22: 272-278
4. Juni P, Dieppe P, Donovan J, Peters T, Eachus J, Pearson N, Greenwood R, Frankel S (2003) Population requirement for primary knee replacement surgery: a cross-sectional study. *Rheumatology (Oxford)* 42: 516-521
5. Merx H, Dreinhöfer KE, Günther KP (2007) Sozialmedizinische Bedeutung der Arthrose in Deutschland. [Socioeconomic relevance of osteoarthritis in Germany]. *Z Orthop Unfall* 145: 421-429
6. Sanders C, Donovan JL, Dieppe PA (2004) Unmet need for joint replacement: a qualitative investigation of barriers to treatment among individuals with severe pain and disability of the hip and knee. *Rheumatology (Oxford)* 43: 353-357
7. Stürmer T, Dreinhöfer K, Gröber-Gratz D, Brenner H, Dieppe P, Puhl W, Günther KP (2005) Differences in the views of orthopaedic surgeons and referring practitioners on the determinants of outcome after total hip replacement. *J Bone Joint Surg Br* 87: 1416-1419

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■ December 10-13
Current Concepts Winter
Orlando, Florida/USA

■ January 31
Symposium de l' Institut
Kerboull 2009
Paris, France

■ February 12-14
ENDO Kongress
Berlin, Germany

■ February 23-24
ORS
Las Vegas, USA

■ February 25-28
76th Annual Meeting American Academy of
Orthopaedic Surgeons (AAOS)
Las Vegas, USA

■ March 13-14
Endoprothetik Forum
Münster, Germany

■ March 20-21
16th Cologne Traumatology Symposium
Cologne, Germany

■ March 20-21
13th AMOU Congress
Halle, Germany

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CeramTec AG