CeraNews

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Torsten Prietzel, MD, Head of the Department for Less-invasive Arthroplasty and Tumour Orthopaedics, Leipzig University Orthopaedic Clinic, Leipzig, Germany; www.minderinvasivendoprothetik.de, mdi-tep@uniklinik-leipzig.de

Prietzel T, Mundt A, Drummer N, Salis-Soglio von G: Less-invasive hip arthroplasty. Report on 850 primary and 34 revision implantations with modified procedure. Poster 58, Annual Meeting of the Orthopaedic Association of Northern Germany, Hamburg, June 18–20, 2009

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Hummel MT, Malkani AL, Yakkanti MR, Baker DL. Decreased Dislocation After Revision Total Hip Arthroplasty Using Larger Femoral Head Size and Posterior Capsular Repair. J Arthroplasty, 24/6, Supplement, 2009

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AAOS
AMERICAN ACADEMY OF
ORTHOPAEDIC SURGEONS

AAOS
AMERICAN ACADEMY OF
ORTHOPAEDIC SURGEONS

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- Shimmin, A. J. et al.:** The effect of operative volume on the outcome of hip resurfacing. Paper No. 316, AAOS Scientific Program 2009
- Aulakh, T. S. et al.:** Learning effect in hip resurfacing: A comparison between developer and trainee surgeons. Paper No. 319, AAOS Scientific Program 2009
- Lavigne, M. et al.:** RCT comparing clinical outcome and gait characteristics after large head THA and hip resurfacing. Paper No. 318, AAOS Scientific Program 2009
- Grammatopoulos, G. A. et al.:** Outcome of THA following revision hip resurfacing is cause dependent and is worse than primary THA. Paper No. 325, AAOS Scientific Program 2009
- Holubowycz, O. et al.:** Large articulation reduces early dislocation after hip replacement: A randomized controlled trial. Paper No. 499, AAOS Scientific Program 2009
- Chavan, P. et al.:** Optimal cup placement in patients with lumbo-sacral deformities. Paper No. 500, AAOS Scientific Program 2009
- Goodman, S. B.:** Clinical performance of modern bearing couples. Symposium on current controversies in bearing surface science, Paper No. II, AAOS Scientific Program 2009
- Walter, W. L.:** Tribology of ceramic bearings: challenges and opportunities. Symposium on current controversies in bearing surface science, Paper No. IV, AAOS Scientific Program 2009
- Berry, D. J. et al.:** Bearing surface selection for total hip arthroplasty: What's best for the young patient? Instructional course 242, AAOS 2009
- Beaule, P. E. et al.:** Technology 2009: Alternative bearing surfaces: The good, bad and indifferent. Instructional course 143, AAOS 2009
- Ranawat, C.:** Ceramic-ceramic THA at 10+ years: What have we learned about what works and what does not? Hip Society/American Association of Hip and Knee Surgeons, AAOS 2009
- Hozack, W.:** Ceramic squeak: how frequent? What causes it? Can it be avoided? Hip Society/American Association of Hip and Knee Surgeons, AAOS 2009
- Learmonth, I.:** Metal-metal bearings in hip arthroplasty: local and systemic biologic responses. Hip Society/American Association of Hip and Knee Surgeons, AAOS 2009

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- ALVAL = aseptic, lymphocytic vasculitis and associated lesions
(無菌性淋巴細胞血管炎與相關病損)
- CoC = ceramic-on-ceramic (陶瓷-陶瓷)
- HRA = hip resurfacing arthroplasty (髖關節表面置換)
- MoPE = metal on PE (金屬-聚乙烯)
- MoM = metal on metal (金屬-金屬)
- PT = pseudo-tumour (假性腫瘤)
- ROM = range of motion (活動範圍)
- THA = total hip arthroplasty (全髖關節置換術)
- UHMWPE = ultra high molecular weight polyethylene
(超高分子量聚乙烯)
- XPE = crosslinked polyethylene (高交聯聚乙烯)

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Instructional Course, Bearing Surface Selection for Total Hip Arthroplasty: What's best for the young patient? Daniel J Berry, MD, Rochester, MN, Steven J MacDonald, MD, London, ON, Canada, Peter F Sharkey, MD, Philadelphia, PA, Jay R Lieberman, MD, Farmington, CT; Papers AAOS 2009

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Hernigou P, Zilber S, Filippini P, Poignard A. Ceramic-Ceramic Bearing Decreases Osteolysis: A 20-year Study versus Ceramic-Polyethylene on the Contralateral Hip. Clin Orthop Rel Res. Published online 13 March 2009

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Harvie P, Giele H, Fang C, Ansoorge O, Ostlere S, Gibbons M, Whitwell D. The

treatment of femoral neuropathy due to pseudotumour caused by metal-on-metal resurfacing arthroplasty. *Hip International* 2008;18 (4): 313–320

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Hart AJ, Skinner JA, Winship P, Faria N, Lukinskaya E, Webster D, Muirhead-Allwood S, Aldam CH, Anwar H, Powell JJ. Circulating levels of cobalt and chromium from metal-on-metal hip replacement are associated with CD8⁺ T-cell lymphopenia. *J Bone Joint Surg Br* 2009; 91-B: 835–42

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Hart AJ, Sabah S, Henckel J, Lewis A, Cobb J, Sampson B, Mitchell A, Skinner JA. The painful metal-on-metal hip resurfacing. *J Bone Joint Surg Br* 2009; 91-B: 738–44

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Hart AJ, Buddhev P, Winship P, Faria N, Powell JJ, Skinner JA. Cup inclination angle of greater than 50 degrees increases whole blood concentrations of cobalt and chromium ions after metal-on-metal hip resurfacing. *Hip International* 2008; 18(3): 212–219

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Peacock A, Say J, Lawrence T. Reactive synovitis following hip resurfacing: a case presentation. *Hip International* 2008; 18(3): 224–7

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- www.freebooks4doctors.com

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- www.freemedicaljournals.com
- www.gfmer.ch/Medical_journals/Free_medical.php
- www.MedicalJournals.co.uk

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- www.freebookcentre.net/medical_text_books_journals/orthopedics_ebooks_online_texts_download.html

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- www.freemedicaljournals.com/fmj/IP_ORTHO.HTM

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- Biomaterials Network: www.biomat.net
- Society for Biomaterials: www.biomaterials.org
- The European Society for Biomaterials: www.esbiomaterials.eu
- Italian Biomaterial Society: www.biomateriali.org
- Material Research Society: www.mrs.org
- Interuniversities Research Center on Materials for Biomedical Engineering: www.cirmib.ing.unitn.it
- The American Ceramic Society: www.ceramics.org

■ **October 11–15**
69th Annual Scientific Meeting (AOA)
Cairns, Australia

■ **October 15–17**
KOA Korean Orthopaedic Association autumn meeting
Seoul, Korea

■ **October 21–24**
DKOU
Berlin, Germany

■ **October 22–24**
ISTA
Hawaii, USA

■ **October 27–30**
39th Annual Advances in Arthroplasty Course
Cambridge, USA

■ **October 29 – November 1**
SICOT
Pattaya, Thailand

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■ **November 24–29**
IOACON – 54th Annual Conference of Indian Orthopaedic Association
Bhubaneswar, India

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CeramTec AG